

KERN FFN15K2IPN



Compact organ balance, protected against water splashes through stainless steel construction and IP65 rating, also with EC type approval [M]



Measuring system

| | |
|---|----------------------|
| Adjustment options: | External calibration |
| Linearity: | 4 g |
| Number of weighing ranges: | 1 |
| Readability [d]: | 2 g |
| Recommended adjusting weight: | 15 kg (M3) |
| Repeatability: | 2 g |
| Resolution: | 7500 |
| Stabilisation time under laboratory conditions: | 2 s |
| Tare range: | -15 kg |
| Warm up time: | 30 min |
| Weighing capacity [Max]: | 15 kg |
| Weighing system: | Strain gauge |
| Weighing units: | kg lb |

Display

| | |
|-------------------------------|----------|
| Display digit height (large): | 2,500 cm |
|-------------------------------|----------|

Environmental conditions

| | |
|--------------------------------|--------|
| Maximum humidity: | 95 % |
| Maximum operating temperature: | 40 °C |
| Minimum ambient temperature: | -10 °C |
| Minimum humidity: | 25 % |

Power supply

| | |
|---------------------------------|------------------------|
| Battery: | no battery |
| Charging time: | 12 h |
| Input voltage: | 110 V - 230 V AC |
| Mains adapter: | Mains adapter external |
| Mains adapter/adapter included: | EURO CH UK US |
| Operating time (Backlight off): | 50 h |
| Operating time (Backlight on): | 30 h |

Service

Adjustment at the location of installation: 961-248

| | |
|--------------------|---------|
| DAkkS Certificate: | 963-128 |
|--------------------|---------|

Category

| | |
|----------------|-------------|
| Brand: | KERN |
| Category: | Balances |
| Product Group: | Bench scale |

Packaging & shipping


| | |
|-------------------------------|--------------------|
| Delivery: | 24 h |
| Dimensions packaging (WxDxH): | 380 x 403 x 206 mm |
| Gross weight: | 4,600 kg |
| Net weight: | 3,850 kg |
| Shipping method: | Parcel service |
| UN number: | 2800 |

Construction

| | |
|-------------------------------------|--------------------|
| Dimensions housing (WxDxH): | 243 x 300 x 127 mm |
| Level indicator: | yes |
| Material weighing plate: | stainless steel |
| Overall dimensions mounted (WxDxH): | 243 x 300 x 127 mm |
| Revolving screw feet: | yes |
| Weighing surface (WxD): | 233 x 193 mm |

Piktogramme

STANDARD



OPTION



FACTORY

