

KERN 313-04

Set of weights, cylindrical, polished stainless steel



Construction

Material:	polished stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	1 mg 200000 mg
OIML class:	E2

Accessories

Wooden box:	313-040-100
-------------	-------------

Service

DAkkS Certificate:	962-303
DAkkS Certificate (Recalibration):	962-303R
Positive tolerance:	921-303
Verification:	952-303

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	150 x 140 x 65 mm
Gross weight:	1,270 kg
Net weight:	1,270 kg
Shipping method:	Parcel service