

KERN 313-05

Set of weights, cylindrical, polished stainless steel



Construction

Material:	polished stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	1 mg 500000 mg
OIML class:	E2

Accessories

Wooden box:	313-050-100
-------------	-------------

Service

DAkkS Certificate:	962-304
DAkkS Certificate (Recalibration):	962-304R
Positive tolerance:	921-304
Verification:	952-304

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	230 x 180 x 110 mm
Gross weight:	3,190 kg
Net weight:	3,190 kg
Shipping method:	Parcel service