

KERN 313-06

Set of weights, cylindrical, polished stainless steel



Construction

Material:	polished stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	1 mg 1000000 mg
OIML class:	E2

Accessories

Wooden box:	313-060-100
-------------	-------------

Service

DAkkS Certificate:	962-305
DAkkS Certificate (Recalibration):	962-305R
Positive tolerance:	921-305
Verification:	952-305

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	230 x 180 x 110 mm
Gross weight:	4,210 kg
Net weight:	4,210 kg
Shipping method:	Parcel service