

KERN 313-07

Set of weights, cylindrical, polished stainless steel



Construction

Material:	polished stainless steel
Shape:	cylindrical

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	1 mg 2000000 mg
OIML class:	E2

Accessories

Wooden box:	313-070-100
-------------	-------------

Service

DAkkS Certificate:	962-306
DAkkS Certificate (Recalibration):	962-306R
Positive tolerance:	921-306
Verification:	952-306

Category

Brand:	KERN
--------	------

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	310 x 180 x 140 mm
Gross weight:	9,010 kg
Net weight:	9,010 kg
Shipping method:	Parcel service