

# KERN 313-08

Set of weights, cylindrical, polished stainless steel



## Construction

|           |                          |
|-----------|--------------------------|
| Material: | polished stainless steel |
| Shape:    | cylindrical              |

## Piktogramme

### STANDARD



### OPTION



## Measuring system

|                |                   |
|----------------|-------------------|
| Nominal value: | 1 mg   5000000 mg |
| OIML class:    | E2                |

## Accessories

|             |             |
|-------------|-------------|
| Wooden box: | 313-080-100 |
|-------------|-------------|

## Service

|                                    |          |
|------------------------------------|----------|
| DAkkS Certificate:                 | 962-307  |
| DAkkS Certificate (Recalibration): | 962-307R |
| Positive tolerance:                | 921-307  |
| Verification:                      | 952-307  |

## Category

|        |      |
|--------|------|
| Brand: | KERN |
|--------|------|

## Packaging & shipping

|                               |                    |
|-------------------------------|--------------------|
| Delivery:                     | 24 h               |
| Dimensions packaging (WxDxH): | 310 x 180 x 140 mm |
| Gross weight:                 | 17,810 kg          |
| Net weight:                   | 17,810 kg          |
| Shipping method:              | Parcel service     |