

# KERN 316-01

Single weights, compact shape, polished stainless steel



## Packaging & shipping

|                               |                |
|-------------------------------|----------------|
| Delivery:                     | 24 h           |
| Dimensions packaging (WxDxH): | 6 x 6 x 6 mm   |
| Gross weight:                 | 1 g            |
| Net weight:                   | 1 g            |
| Shipping method:              | Parcel service |

## Construction

|             |                          |
|-------------|--------------------------|
| Dimensions: | 6   6                    |
| Material:   | polished stainless steel |
| Shape:      | compact shape            |

## Piktogramme

### STANDARD



### OPTION



## Measuring system

|                |          |
|----------------|----------|
| Nominal value: | 1 g      |
| OIML class:    | E2       |
| Tolerance:     | 0,030 mg |

## Accessories

|                |             |
|----------------|-------------|
| Aluminium box: | 317-010-600 |
| Plastic box:   | 317-020-400 |
| Wooden box:    | 317-010-100 |

## Service

|                                    |          |
|------------------------------------|----------|
| DAkKS Certificate:                 | 962-331  |
| DAkKS Certificate (Recalibration): | 962-331R |
| Positive tolerance:                | 921-331  |
| Verification:                      | 952-331  |

## Category

|        |      |
|--------|------|
| Brand: | KERN |
|--------|------|