

KERN 316-02

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	6 x 6 x 10 mm
Gross weight:	2 g
Net weight:	2 g
Shipping method:	Parcel service

Construction

Dimensions:	6 10
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	2 g
OIML class:	E2
Tolerance:	0,040 mg

Accessories

Aluminium box:	317-020-600
Plastic box:	317-020-400
Wooden box:	317-020-100

Service

DAkKS Certificate:	962-332
DAkKS Certificate (Recalibration):	962-332R
Positive tolerance:	921-332
Verification:	952-332

Category

Brand:	KERN
--------	------