

KERN 316-03

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	8 x 8 x 13 mm
Gross weight:	10 g
Net weight:	10 g
Shipping method:	Parcel service

Construction

Dimensions:	8 13,200
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	5 g
OIML class:	E2
Tolerance:	0,050 mg

Accessories

Aluminium box:	317-030-600
Plastic box:	317-030-400
Wooden box:	317-030-100

Service

DAkKS Certificate:	962-333
DAkKS Certificate (Recalibration):	962-333R
Positive tolerance:	921-333
Verification:	952-333

Category

Brand:	KERN
--------	------