

KERN 316-04

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	10 x 10 x 17 mm
Gross weight:	10 g
Net weight:	10 g
Shipping method:	Parcel service

Construction

Dimensions:	10 17,300
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	10 g
OIML class:	E2
Tolerance:	0,060 mg

Accessories

Aluminium box:	317-040-600
Plastic box:	317-040-400
Wooden box:	317-040-100

Service

DAkKS Certificate:	962-334
DAkKS Certificate (Recalibration):	962-334R
Positive tolerance:	921-334
Verification:	952-334

Category

Brand:	KERN
--------	------