

KERN 316-05

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	13 x 13 x 20 mm
Gross weight:	20 g
Net weight:	20 g
Shipping method:	Parcel service

Construction

Dimensions:	13 20
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	20 g
OIML class:	E2
Tolerance:	0,080 mg

Accessories

Aluminium box:	317-050-600
Plastic box:	317-050-400
Wooden box:	317-050-100

Service

DAkKS Certificate:	962-335
DAkKS Certificate (Recalibration):	962-335R
Positive tolerance:	921-335
Verification:	952-335

Category

Brand:	KERN
--------	------