

KERN 316-06

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	18 x 18 x 26 mm
Gross weight:	50 g
Net weight:	50 g
Shipping method:	Parcel service

Construction

Dimensions:	18 26,200
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	50 g
OIML class:	E2
Tolerance:	0,100 mg

Accessories

Aluminium box:	317-060-600
Plastic box:	317-060-400
Wooden box:	317-060-100

Service

DAkKS Certificate:	962-336
DAkKS Certificate (Recalibration):	962-336R
Positive tolerance:	921-336
Verification:	952-336

Category

Brand:	KERN
--------	------