

KERN 316-08

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	28 x 28 x 43 mm
Gross weight:	200 g
Net weight:	200 g
Shipping method:	Parcel service

Construction

Dimensions:	28 43
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	200 g
OIML class:	E2
Tolerance:	0,300 mg

Accessories

Aluminium box:	317-080-600
Plastic box:	317-080-400
Wooden box:	317-080-100

Service

DAkKS Certificate:	962-338
DAkKS Certificate (Recalibration):	962-338R
Positive tolerance:	921-338
Verification:	952-338

Category

Brand:	KERN
--------	------