

KERN 316-09

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	38 x 38 x 58 mm
Gross weight:	500 g
Net weight:	500 g
Shipping method:	Parcel service

Construction

Dimensions:	38 58,500
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	500 g
OIML class:	E2
Tolerance:	0,800 mg

Accessories

Aluminium box:	317-090-600
Plastic box:	317-090-400
Wooden box:	317-090-100

Service

DAkKS Certificate:	962-339
DAkKS Certificate (Recalibration):	962-339R
Positive tolerance:	921-339
Verification:	952-339

Category

Brand:	KERN
--------	------