

KERN 316-11

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	48 x 48 x 73 mm
Gross weight:	1 kg
Net weight:	1 kg
Shipping method:	Parcel service

Construction

Dimensions:	48 73
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	1e+006 mg
OIML class:	E2
Tolerance:	1,600 mg

Accessories

Aluminium box:	317-110-600
Plastic box:	317-110-400
Wooden box:	317-110-100

Service

DAkKS Certificate:	962-341
DAkKS Certificate (Recalibration):	962-341R
Positive tolerance:	921-341
Verification:	952-341

Category

Brand:	KERN
--------	------