

KERN 316-12

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	60 x 60 x 94 mm
Gross weight:	2 kg
Net weight:	2 kg
Shipping method:	Parcel service

Construction

Dimensions:	60 94
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	2e+006 mg
OIML class:	E2
Tolerance:	3 mg

Accessories

Aluminium box:	317-120-600
Plastic box:	317-120-400
Wooden box:	317-120-100

Service

DAkkS Certificate:	962-342
DAkkS Certificate (Recalibration):	962-342R
Positive tolerance:	921-342
Verification:	952-342

Category

Brand:	KERN
--------	------