

KERN 316-13

Single weights, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	80 x 80 x 132 mm
Gross weight:	5 kg
Net weight:	5 kg
Shipping method:	Parcel service

Construction

Dimensions:	80 132,500
Material:	polished stainless steel
Shape:	compact shape

Piktogramme

STANDARD



OPTION



Measuring system

Nominal value:	5e+006 mg
OIML class:	E2
Tolerance:	8 mg

Accessories

Aluminium box:	317-130-600
Plastic box:	317-130-400
Wooden box:	317-130-100

Service

DAkkS Certificate:	962-343
DAkkS Certificate (Recalibration):	962-343R
Positive tolerance:	921-343
Verification:	952-343

Category

Brand:	KERN
--------	------